## Gebrauchsanweisung Implantate

### Implants instruction for use

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## Implants instruction for use

**Experience in surgery and basic knowledge in implantology are prerequisites for using the implant system. In any case we do recommend to follow broader theoretical and practical knowledge in the field of implantology through studies and workshops. The sole knowledge of this instruction for use is not sufficient for the usage of this system!**

#### GENERELL PRECAUTIONS

- Dentegris-Implants may only be used by qualified medical professionals with respectvely applicable training.
- Components that are not part of the system produced by Dentegris, may have a negative effect on the function and effectiveness of the implants.
- Non sterile and used implants must not be used again. Dentegris is not liable for the usage of such implants even if sterilised by the medical professional and does not accept the return of implants if the seal of the packaging is broken.
- The implants must be stored in their original packaging in a dry and from sunlight sheltered location. Before the implant is used it has to be checked if the packaging and the implant are undamaged.
- The implant must not be used if the expiry date of the sterilisation has been exceeded.
- Information on the usage of Dentegris-Implants is available in the surgical protocol and the product catalog, which contain the current knowledge of Dentegris. This information is a recommendation only and cannot substitute the duty of the user to assure that the product is suitable for the intended use.
- The usage and application of the products is beyond our control and Dentegris is not liable for any damages that may result thereof.
- Dentegris is not liable, if the above instructions are not obeyed.

#### AREAS OF USE

- The drills, countersink and/or tab of the respective implant diameter must be used for the preparation of the insertion of the implant.
- Drills as well as countersink and tab are to be applied with even force and sufficient sterile cooling.
- The planing of the implantation should be supported by applicable x-ray diagnostics. Complicated cases should rely on the latest 3D diagnostics in order to avoid failures.
- Dentegris mini implants should only be used for stabilizing removable full dentures in the lower jaw if the sagittal width of the jaw bone is insufficient for Dentegris implants with diameters from 3.3 mm. At least four parallel and equidistant Dentegris mini implants have to be placed in the interforaminal area with a minimum space of 5 mm between the implant sites to stabilize a full denture of the lower jaw. Any inclination of the individual implants toward one another should not exceed 8 degrees.
- Dentegris implants are suited for endosseous placement in the upper and lower jaw with one-phase or two-phase surgical techniques. All indication areas for endosseous oral implants in the upper and lower jaw (treatment of single tooth, bridge abutments, partial full dentures, orthodontics) apply. If there is sufficient bone supply, it is recommended to use implants with diameters 4.1, 4.5 and 5.5; if the bone supply is less, diameters 3.3 and 3.8 are better suited. The following applies to all indications: the diameter and length of the implants must be in a balanced proportion to the clinical crown. Implants should be centered underneath the prosthetic crown. (No axial load in the outlying area of the implant platform).

Up- per jaw		3,3	-	(+)	-	-	-	-	-	-
		3,8	+	++	+	+	+	+	(+)	(+)
		4,1	++	-	++	++	++	+	+	+
		4,5	++	-	++	++	++	++	++	++
		5,5	++	-	++	++	++	++	++	++
	<b>Tooth</b>	1	2	3	4	5	6	7		
Lo- wer jaw		3,3	(+)	(+)	-	-	-	-	-	-
		3,8	++	++	++	+	+	+	(+)	(+)
		4,1	-	-	++	++	++	+	+	+
		4,5	-	-	++	++	++	++	++	++
		5,5	-	-	++	++	++	++	++	++

-	= not suitable
(+)	= suitable, in conjunction with teeth or other implants only
+	= suitable, preferably in conjunction with teeth or other implants
++	= suitable, including for individual crown without interconnection

- In the case of mechanical insertion of the implant, this should not exceed 35 Ncm for the Dentegris-Mini-Implants and 50 Ncm for the Dentegris-Implants.
- For taking the impression, the respective impression posts and analogs of the corresponding Dentegris system are to be used.
- For the prosthetics of the Dentgris Implant System and the Dentegris Mini-Implant System, the respective abutments of the corresponding system are to be used.

#### CONTARINDICATIONS

The scientific data on absolute and relative contraindications in implantology have to be taken into account.

#### SIDE EFFECTS

In order to reduce side effects to a minimum, the patient should not physically exercise after the operation. The area of the insertion may be subject to local pain and swolleness after the operation. In rare cases chronic pain can remain as a result of an implantation as well as reduction of the bone in the area of the alveolar in the long run. Very rare complications are fractures of the jaw and damages to nerves, which can be prevented by careful preoperative diagnostics and atraumatic surgical procedure. Dentegris does recommend informing the patient on the risks of implantology and the surgical procedures necessary in writing as well as through a personal consultation.

#### PACKAGING AND STERILITY

Dentegris-Implants are delivered in a double-sterile blister, which guarantees its sterility, if the packaging is undamaged and correctly stored in a protected location. The packaging composes of the following elements:

- The carton with instructions for use and sticker for the patient's file
- Double-sterile- blister with labels and seals to guarantee sterility of the implant; (inside are the implant and carrier with cover screw).

#### HOW TO TAKE THE IMPLANT OUT OF THE STERILE PACKAGING

The regulations for working under sterile conditions must be followed when removing the implant from the packaging. The packaging must not be opened prior to the surgical insertion, but when the implant bed has been fully prepared. In any case the contact with non-sterile items must be avoided.

#### REMOVE THE IMPLANT FROM ITS PACKAGING AS FOLLOWS:

- Assistant personal opens the outer blister and sterile assistant personal takes out the inner blister from it.
- Sterile personal opens the inner blister.
- The implantologist takes out the implant and inserts it immediately into the prepared part of the jaw.
- When sufficient friction has been has been reached remove the handle and unscrew the cover screw from it.
- Insert the implant into its final position with the chosen mechanical instrument.

#### DOCUMENTATION AND IDENTIFICATION

Dentegris recommends a systematic pre- and post-operational x-ray documentation as well as a careful clinical documentation. The implants can be identified through the article code and the lot number on the packaging and these identification numbers must be documented in the patients file or surgical report.

**If the packaging or part of it is damaged or the sterility has expired that implant must not be used under any circumstances.**

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