MucoMatrixX®

DE	GEBRAUCHSINFORMATION
EN	INSTRUCTIONS FOR USE
FR	INSTRUCTIONS D'UTILISATION
ES	INSTRUCCIONES DE USO
IT	ISTRUZIONI PER L'USO
PT	INFORMAÇÕES DE UTILIZAÇÃO
HU	HASZNÁLATI ÚTMUTATÓ
HR	INFORMACIJA O KORIŠTENJU
TR	KULLANIM BILGILERI
RU	ИНФОРМАЦИЯ ПО ПРИМЕНЕНИЮ
EL	ΠΛΗΡΟΦΟΡΙΕΣ ΧΡΗΣΗΣ

CE 0483

Fabricant / Manufacturer
MBP
Medical Biomaterial
Products GmbH
Lederstraße 7
D-19306 Neustadt-Glewe

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Instructions for use MucoMatrixX®**MucoMatrixX® composition and pack size**

- Porcine collagen, manufactured from porcine corium
- Single use
- Thickness varies between 0.8 and 5.0 mm
- MucoMatrixX® is offered in a sterile double packaging

The following sizes are available:

- 15 x 20 mm
- 20 x 30 mm
- 30 x 40 mm

Therapeutic indications

MucoMatrixX® is indicated for dental regeneration of soft tissue defects and regenerating soft tissues around teeth or implants.

- Treatment of recession defects combined with all the current mucogingival surgical techniques (e.g. coronal and lateral advancement flap, tunnel technique)
- Soft tissue augmentation around teeth and implants to widen the keratinised gingiva (vestibuloplasty)
- For surgical bone and bone wall defects or in the scope of jaw ridge augmentation combined with bone replacement materials where simultaneous thickening of the soft tissue is desired
- Covering implants after immediate or delayed implantation for extraction alveoli
- Jaw ridge regeneration for prosthetic restorations
- Closing extraction alveoli for socket preservation (socket seal technique)
- Soft tissue regeneration of the jaw ridge in preparation for subsequent augmentative interventions

Contraindications

MucoMatrixX® must not be used in patients with

- Acute infection of the buccal cavity, or acute or chronic inflammation at the implantation site

- General illness which do not allow for oral and maxillofacial surgery, implantological, periodontal, endodontic or other oral surgical measures
- Known hypersensitivity to porcine collagen

Use during pregnancy and breast-feeding

There is no information available on the use of MucoMatrixX® during pregnancy and breastfeeding, or its impact on human fertility. The treating dentist must therefore make an individual assessment of the benefit to the mother and the possible risks for the child before using MucoMatrixX®.

Use in children and elderly patients

There is no information available indicating special precautions necessary pertaining to the age of the patient.

Instructions for use

MucoMatrixX® is only to be used by dentists familiar with controlled bone and tissue regeneration techniques based on the appropriate training. MucoMatrixX® is of natural origin and therefore slight undulation of the collagen structure and variations in the membrane thickness may occur in its dry state. These appearances do not impact the quality and functionality of MucoMatrixX®.

While maintaining its sterility, MucoMatrixX® may be cut to the shape and size of the defect to be treated using scissors or a scalpel, preferably in a rehydrated state. Using suitable templates may be helpful when determining the area of MucoMatrixX® required. The matrix may be affixed by sutures or pins and may be recommended depending on the application to prevent dislocation or micro-movements, thus promoting faster healing. MucoMatrixX® may be used for both closed healing wounds and for open healing. During open healing, the supply to and revascularisation of the matrix must be guaranteed, e.g. through close contact with the underlying periosteum. Always avoid exposure when used in recession coverage, (see 'Special instructions for use in treating gingival recession defects').

The basic principles of aseptic technique and patient medication must be observed when using MucoMatrixX®.

- The necessary surgical intervention begins with exposing the defect
- Where clinically indicated, the bone defect is filled with a suitable augmentation material (for example autogenous bone, allogenic, xenogenic or alloplastic bone replacement materials)
- The outer packaging, which is sterile inside, should be removed, for example by an assistant in the unsterile surgical area, and opened under sterile conditions. The inner packaging, which is sterile on both the outside and inside, should be removed under sterile conditions
- Only apply and trim MucoMatrixX® when moist. For this purpose the membrane must be rehydrated for approximately 5–20 minutes in sterile saline solution or blood before application, depending on the type of use and the desired flexibility of the matrix
- The shape and size of the matrix should be adapted to the size of the defect. Following rehydration, MucoMatrixX® may be cut to the desired shape by scissors or a scalpel. With reduced rehydration times rounding off the corners may help to prevent perforation of the gingival soft tissue after flap closure. MucoMatrixX® is placed over the defect and adapted in situ using moderate pressure. The amount of time during which pressure is applied varies by the extent of bleeding. Adhesion to the bone surface occurs through swelling and gel formation of the collagen fibres brought about by contact with blood
- Completely soaking MucoMatrixX® in blood and exudates allows for anchoring to the defective structure or augmentation material, as well as forming a blood clot under the matrix
- When preparing a split flap, the MucoMatrixX® must be sutured to the intact periosteum to ensure close contact between the matrix and the periosteal wound bed, and immobilisation must be ensured. Single button sutures or cruciate sutures may be used; the use of resorbable sutures is recommended
- Following surgical intervention the patient should follow dental hygiene per the dentist's instructions

Special instructions for treatment of gingival recession defects

On principle, predictable coverage of the recession may only be expected for treatment of Miller Class I and II defects.

Prolonged rehydration of 10–20 minutes is recommended in order to guarantee sufficient flexibility of the matrix when applying

MucoMatrixX® using the tunnel technique. The matrix is also recommended to be affixed with simple interrupted sutures or cruciate sutures when using the tunnel technique. Exposure should always be avoided when using MucoMatrixX® for gingival recessions. Ensure the repositioned flaps fully cover the matrix. Primary wound closure allows blood vessels to penetrate and incorporate into the matrix. Premature exposure may lead to transplant loss. Avoid any form of mechanical trauma to the area treated after the operation. The patient should refrain from brushing teeth in the affected area for 4 weeks. Plaque may be prevented using a mouthwash with 0.12 % chlorhexidine solution. The patient should be seen weekly following the operation to assess recovery and to check the formation of plaque.

Adverse effects

- In rare cases, allergic reactions to collagen membranes cannot be excluded
- In very rare cases, intolerance to collagen may occur
- In rare cases the prolonged resorption time may lead to inflammatory tissue reactions
- As with any foreign material, existing local infections may be worsened by implanting MucoMatrixX®
- Possible general complications may result from the surgical intervention itself, for example, gingival recession, severe gum bleeding, swelling of the soft tissue, temperature sensitivity, desquamation of the gingiva epithelium in the flap area

Interactions with other products

- Do not use antiseptics releasing chlorine (e.g. chloramine), as well as tannin and caustics which modify proteins in combination with collagen
- Ointments and powders as well as silicon preparations may alter the spaces between collagen fibres, fleeces or sponges, so these should also not be used concomitantly

Warnings and precautions

The resorption time may be accelerated in open healing or MucoMatrixX® exposure during the healing phase.

MucoMatrixX® is only indicated for the specified applications. The attending dentist must educate patients about possible contraindications, adverse effects and necessary precautions. In the event of post-operative complaints, for example pain, infection or other unusual symptoms, the patient should immediately consult a dentist. Patients with serious general medical conditions (e.g. poorly controlled diabetes mellitus, serious blood pressure issues, severe peripheral artery occlusive disease (PAOD), malignant tumours or autoimmune diseases) or patients undergoing for example long-term steroid therapy or anti-coagulant therapy, must be treated with special care with respect to all surgical interventions. Do not reuse MucoMatrixX® once removed from the packaging and/or following patient contact, as this presents an increased risk of contamination with subsequent risk of infection.

Shelf life

The expiry date is imprinted on the box and sterile inner packaging. MucoMatrixX® must not be used after the expiration date.

Disposal

After use, handle and dispose of all the unused product and packaging according to recognised medical practices and the applicable national and regional environmental protection laws for disposal of packaging materials and biological materials.

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